

Streamer Feedback Validation (Game/Daniil's Tasks)

Vicon and Streamer Setup (1)

- Connect Laptop to TV
- Start and Setup the game
- Wait for Vicon Setup
- Start Streamer
- Select Faceplate Device
- Select Hip Angles
- Enter Leg Length
- Select Moving Average Filter (Size 5)
- Ask to stand straight and then Zero
- Select Device to Plot
- Enter Subject Name
- Enter save Path
- Enter IP Address

Recording Sequence (2, 4-15)

➤ Start Screen Recording	■												■
➤ Start Vicon Recording													
➤ Start Streamer Recording													
➤ Start Treadmill													
➤ Stop Treadmill													
➤ Stop Streamer Recording													
➤ Stop Vicon Recording													
➤ Stop Screen Recording	■												■

Targets and Worm Setup (3)

- Process First Trial and get Targets
- Write down targets:
 - a. Slider Max: _____
 - b. Slider Target: _____ ± _____
 - c. Slider Min: _____
- Start Streaming
- Enter Worm Parameters (NEGATED)
- Ask to Walk in Place and Check Worm is moving
- Turn on TV